



User's Manual

PI3-kinase p85-alpha/gamma ELISA Kit



DEIA-XYA1439



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PI3-Kinase p85-alpha/gamma Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PI3-Kinase p85-alpha/gamma and the effects that certain stimulation conditions have on PI3-Kinase p85-alpha/gamma expression in different cell lines.

Storage

4°C/6 Months