



User's Manual

60S Ribosomal Protein L10 ELISA Kit



DEIA-XYA144



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The 60S Ribosomal Protein L10 Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can monitor 60S Ribosomal Protein L10 protein expression profile in cells. The kit can be used for measuring the relative amounts of 60S Ribosomal Protein L10 in cultured cells as well as screening for the effects that various treatments, inhibitors (ie. siRNA or chemicals), or activators have on 60S Ribosomal Protein L10.

Storage

4°C/6 Months