



## User's Manual

# PKC alpha ELISA Kit



DEIA-XYA1451



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The PKC alpha Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PKC alpha and the effects that certain stimulation conditions have on PKC alpha expression in different cell lines.

### Storage

4°C/6 Months