



User's Manual

PKD1/PKC mu ELISA Kit

REF

DEIA-XYA1470



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PKD1/PKC mu Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PKD1/PKC mu and the effects that certain stimulation conditions have on PKD1/PKC mu expression in different cell lines.

Storage

4°C/6 Months