



User's Manual

PLCbeta3 ELISA Kit

REF

DEIA-XYA1478



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PLCbeta3 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PLCbeta3 and the effects that certain stimulation conditions have on PLCbeta3 expression in different cell lines.

Storage

4°C/6 Months