



User's Manual

POU4F1 ELISA Kit

REF

DEIA-XYA1496



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The POU4F1 DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor POU4F1 in eukaryotic nuclear or cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The POU4F1 DNA-Binding ELISA Kit detects endogenous levels of POU4F1.