



User's Manual

PPAR-alpha ELISA Kit

REF

DEIA-XYA1499



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The PPAR-alpha DNA-Binding ELISA Kit contains components necessary for qualitative detection of active transcription factor PPAR-alpha in eukaryotic nuclear or cell lysates.

Storage

4°C/6 Months

Reference Values



Specificity

The PPAR-alpha DNA-Binding ELISA Kit detects endogenous levels of total PPAR-alpha protein.