



User's Manual

PRKAG1/2/3 ELISA Kit

REF

DEIA-XYA1506



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PRKAG1/2/3 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of PRKAG1/2/3 and the effects that certain stimulation conditions have on PRKAG1/2/3 expression in different cell lines.

Storage

4°C/6 Months