



User's Manual

Progesterone Receptor ELISA Kit



DEIA-XYA1507



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Progesterone receptor Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Progesterone receptor and the effects that certain stimulation conditions have on Progesterone receptor expression in different cell lines.

Storage

4°C/6 Months