



User's Manual

Prostate Apoptosis Response protein-4 ELISA Kit



DEIA-XYA1512



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Prostate Apoptosis Response protein-4 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Prostate Apoptosis Response protein-4 and the effects that certain stimulation conditions have on Prostate Apoptosis Response protein-4 expression in different cell lines.

Storage

4°C/6 Months



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com