



## User's Manual

# Proteasome alpha3 ELISA Kit



DEIA-XYA1513



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Proteasome alpha3 Colorimetric Cell-Based ELISA Kit allows for the qualitative detection of Proteasome alpha3 and the effects that certain stimulation conditions have on Proteasome alpha3 expression in different cell lines.

### Storage

4°C/6 Months