



User's Manual

Rat GRO/KC (Colorimetric) ELISA Kit



DEIA-XYA1531



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The GRO/KC (Rat) contains the components necessary for quantitative determination of natural or recombinant Rat GRO/KC concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

16-1000 pg/mL

Specificity

The Rat GRO/KC ELISA is capable of recognizing both recombinant and naturally produced Rat GRO/KC proteins. The antigens listed below were tested at 50 ng/mL and exhibited 100% cross reactivity.

Mouse: KC

The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity or interference.

Human: GRO/MGSA, GRObeta, GROgamma

Mouse: MIP-2

Rat: GRObeta/MIP-2