



User's Manual

Rat MIP-1alpha (Colorimetric) ELISA Kit



DEIA-XYA1539



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The MIP-1alpha (Rat) contains the components necessary for quantitative determination of natural or recombinant Rat MIP-1alpha concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Detection Range

16-1000 pg/mL

Specificity

The Rat MIP-1alpha ELISA is capable of recognizing both recombinant and naturally produced Rat MIP-1alpha proteins. The antigens listed below were tested at 50 ng/mL and exhibited less than 1% cross reactivity.

Mouse: MIP-1alpha

The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: MIP-1alpha, MIP-1beta, IP-10, LAG-1, LD-78

Mouse: MIP-1beta, MIP-1gamma

Rat: MIP-1beta, MCP-1, RANTES