



## User's Manual

# Retinoic Acid Receptor beta ELISA Kit

REF

DEIA-XYA1568



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Retinoic Acid Receptor  $\beta$  Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can monitor Retinoic Acid Receptor  $\beta$  protein expression profile in cells. The kit can be used for measuring the relative amounts of Retinoic Acid Receptor  $\beta$  in cultured cells as well as screening for the effects that various treatments, inhibitors (ie. siRNA or chemicals), or activators have on Retinoic Acid Receptor  $\beta$ .

### Storage

4°C/6 Months