



## User's Manual

# G-CSFR ELISA Kit



DEIA-XYA840



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION


### Intended Use

The Human G-CSFR ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human G-CSFR concentrations within any experimental sample including cell lysates, serum and plasma.

### Storage

4°C/6 Months

### Typical Standard Curve



### Detection Range

32-2000 pg/mL

### Sensitivity

The Human G-CSFR ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human G-CSFR proteins within the range of 32-2000 pg/mL.

### Specificity

The Human G-CSFR ELISA is capable of recognizing both recombinant and naturally produced Human G-

CSFR proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross reactivity or interference.

Human: G-CSF, GM-CSF R $\beta$ , M-CSF R/Fc Chimera

Mouse: G-CSF