



## User's Manual

# Calsenilin/KCNIP3 ELISA Kit

REF

DEIA-XYA873



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Calsenilin/KCNIP3 Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can monitor Calsenilin/KCNIP3 protein expression profile in cells. The kit can be used for measuring the relative amounts of Calsenilin/KCNIP3 in cultured cells as well as screening for the effects that various treatments, inhibitors (ie. siRNA or chemicals), or activators have on Calsenilin/KCNIP3.

### Storage

4°C/6 Months



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)