



## User's Manual

# Apomorphine ELISA Kit

REF

DEIA-XYL7



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

CD's Apomorphine/Apocodeine ELISA (Enzyme-Linked ImmunoSorbent Assay) test kit is a qualitative one-step kit designed for use as a screening device for the detection of drugs and/or their metabolites. The kit was designed for screening purposes and is intended for forensic use only.

### Sensitivity

Apomorphine: 38.3 ng/mL

Apocodeine: 3.9 ng/mL

### Specificity

Apocodeine: 825%

R(-)Apomorphine: 100%

R(-)Norapomorphine: 17.5%

S(+)Apomorphine: 9.54%

R(-)Propylnorapomorphine: 0.88%

Boldine: 0.34%

Imipramine: 0.15%

(+)Isocorydine: 0.14%

Codeine: 0.05%

Levallorphan: 0.03%

Hydromorphone: 0.02%