



User's Manual

Ofloxacin ELISA Kit



DEIA-XYZ34



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Test Kit enables the detection of ofloxacin as low as 0.1 ppb level in meat and tissue samples. This kit provides an accurate, rapid, cost effective and easy to use extraction protocol.

General Description

Ofloxacin is a second-generation fluoroquinolone that inhibits the activity of bacterial DNA gyrase.

Storage

2-8°C

Reconstitution And Storage

The shelf life is 12 months when the kit is properly stored.

Detection Range

The ELISA Test Kit has the capacity for 96 determinations or testing of 42 samples in duplicate.