



**User's Manual**

# **Brombuterol competitive colorimetric ELISA Kit**



**DEIA-XYZ37**



**96T**





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This ELISA Kit offers rapid, quantitative analysis of brombuterol in feed, liver, kidney, milk, serum, plasma and urine samples. This kit is based on a competitive colorimetric ELISA assay.

### Storage

2-8°C

### Reconstitution And Storage

The shelf life is 12 months when the kit is properly stored.

### Detection Range

This ELISA Test Kit has the capacity for 96 determinations or testing of 42 samples in duplicate (assuming 12 wells for standards).

### Detection Limit

Meat & organs: 0.1 ppb

Urine: 0.2 ppb

Feed: 1 ppb