



User's Manual

Dimetridazole ELISA Kit

REF

DEIA-XYZ41



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Test Kit is a competitive colorimetric ELISA assay that is accurate, rapid, cost effective and easy to use and offers a unique extraction method. The Kit is an effective screening tool for detecting dimetridazole in various sample types at concentrations as low as 0.5 ppb.

Storage

2-8°C

Reconstitution And Storage

The shelf life is 12 months when the kit is properly stored.

Detection Range

The kit has the capacity for 96 determinations or testing of 42 samples in duplicate

