



User's Manual

Doramectin ELISA Kit



DEIA-XYZ42



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Test Kit for bovine muscle employs a unique extraction method, which makes the kit accurate, rapid, cost effective and easy to use. The Kit for bovine muscle is a competitive colorimetric ELISA assay.

Storage

2-8°C

Reconstitution And Storage

The shelf life is 12 months when the kit is properly stored.

Detection Range

The kit has the capacity for 96 determinations or testing of 42 samples in duplicate

Detection Limit

5.0 ppb

