



User's Manual

Human HBeAg ELISA Kit



DEIA003



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human HBeAg ELISA Kit is for qualitative detection of hepatitis B virus e antigen (HBeAg) in human serum or plasma. It is intended for use in clinical laboratories for diagnosis and management of patients related to infection with hepatitis B virus.

Storage

Store unopened kit and all kit components at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

Human Anti-HBeAg ELISA Kit