



User's Manual

Human Anti-HBcAg IgM ELISA Kit



DEIA006



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human Anti-HBcAg IgM ELISA Kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative determination of IgM class antibodies to hepatitis B virus core antigen in human serum or plasma. It is intended for use in clinical laboratories for management of patients related to infection with hepatitis B virus.

Storage

The components of the kit will remain stable through the expiration date indicated on the label and package when stored at 2-8°C, do not freeze. For more detailed information, please download the following document on our website.

Sensitivity

Human Anti-HBcAb IgM ELISA Kit