



User's Manual

Anti-human-tissue- transglutaminase IgA ELISA Kit

REF

DEIA008J



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The described assay is intended for the quantitative determination of underglycosylated Bone Sialoprotein in serum.

Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

Precision

Intra-assay Precision: 5.2%

Inter-assay Precision: 13.7%-15.4%