



**User's Manual**

# **Human Anti-Toxoplasma Gondii IgM Antibody, Anti- TOXO IgM ELISA Kit**

**REF**

**DEIA010**



**96T**

**RUO**

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This assay is an enzyme-linked immunosorbent assay (ELISA) designed for the presumptive qualitative detection of IgM antibody to TOX in human serum and for the presumptive diagnosis of acute, recent, or reactive TOX infection.

### Storage

Store the kit at 2-8°C. Keep microwells sealed in a dry bag with desiccants. For more detailed information, please download the following document on our website.

### Sensitivity

The results yielded a Sensitivity of 89.5%.