



User's Manual

C-reactive Protein ELISA Kit

REF DEIA012J

Σ 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The assay is a sandwich enzyme immuno assay intended for the quantitative determination of C-reactive Protein in dried blood samples.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-Assay Precision: 5.5%

Inter-Assay Precision: 11.6%-13.8%

Sensitivity

0.921 ng/mL