



User's Manual

Osteoprotegerin ELISA Kit



DEIA015J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The assay is a sandwich ELISA intended for the determination of mouse/rat Osteoprotegerin in serum, plasma, urine and cell culture supernatants.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com