



User's Manual

Chloramphenicol, CAP ELISA Kit

REF

DEIA018



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Chloramphenicol ELISA Kit is a competitive enzyme immunoassay for the quantitative analysis of Chloramphenicol honey/tissue, milk/serum and feed.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0.025-6.4 ppb

Detection Limit

Honey, Tissue: 0.1 ppb

Milk, Serum: 0.025 ppb

Feed: 100 ppb

Sensitivity

0.025 ppb