



User's Manual

TNF α -Blocker ADA, Antibodies against etanercept ELISA Kit

REF

DEIA018J



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The here described Enzyme-Linked-Immuno-Sorbent-assay (ELISA) Kit is intended for the qualitative determination of free human antibodies against therapeutic anti-TNF α antibodies in the rheumatism therapy in EDTA-plasma and serum.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.