



User's Manual

TNF α -Blocker ADA, Total Antibodies against infliximab ELISA Kit

REF

DEIA020J



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Enzyme Linked Immuno Sorbent assay (ELISA) is intended for the determination of human antibodies against TNF α blocker Infliximab in the presence of Infliximab in EDTA plasma and serum.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-Assay Precision: 3.18%-4.67%

Inter-Assay Precision: 11.41%-12.61%