



User's Manual

Tumor Necrosis Factor alpha [TNF α] ELISA Kit



DEIA027J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The assay is intended for the quantitative determination of TNF- α in plasma, serum, stool and cell culture supernatant.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-Assay Precision: 6.3%

Inter-Assay Precision: 8.2%

Sensitivity

10 pg/mL

Specificity

GM-CSF < 10%

M-CSF < 10%

G-CSF < 10%