



User's Manual

Vitamin D binding Protein (VDBP) ELISA Kit



DEIA028J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay is an enzyme immunoassay intended for the quantitative determination of free and not actin complex bound vitamin D binding protein in serum, plasma and urine.

Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

Precision

Intra-assay Precision: 3.2%-5.0%

Inter-assay Precision: 12.7%

Sensitivity

1.23 ng/mL