



User's Manual

Melamine ELISA Kit



DEIA030



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Melamine ELISA Test Kit is a competitive enzyme immunoassay for the quantitative analysis of Melamine in milk, milk powder, infant formula, meat, fish, feed, chocolate and other food samples.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

10-810 ppb

Detection Limit

Liquid milk (method 1), eggs: 100 ppb

Liquid milk (method2): 80 ppb

Chicken: 200 ppb

Milk, feed: 300 ppb

Sensitivity

10 ppb