



User's Manual

Sulfadimethoxine ELISA Kit



DEIA034



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Sulfadimethoxine ELISA Test Kit is a competitive enzyme immunoassay for the quantitative analysis of Sulfadimethoxine in tissue, egg, honey, serum, urine, milk.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

1-81 ppb

Detection Limit

tissue, egg ,honey: 1 ppb

serum, urine: 4 ppb

Milk: 20 ppb

Sensitivity

1 ppb