



User's Manual

Colistin ELISA Kit



DEIA042



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Colistin ELISA Kit is a competitive enzyme immunoassay for the quantitative analysis of colistin in meat (beef, pork and chicken) samples.

Storage

Reagents must be stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

Meat (Beef/Pork/Chicken): 4 ppb

Egg: 50 ppb