



User's Manual

Sarafloxacin ELISA Kit



DEIA044



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Sarafloxacin ELISA Test Kit is a competitive enzyme immunoassay for the quantitative analysis of Sarafloxacin in aquatic product.

Storage

Reagents must be stored at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0.2-16.2 ppb

Sensitivity

0.2 ppb



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com