



User's Manual

Scl-70 IgG ELISA Kit



DEIA04950



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Scl-70 ELISA Kit is intended for the detection of IgG antibody to Scl-70 in human serum or plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.