



User's Manual

Hemoglobin/Haptoglobin- Complex (1-point calibration) ELISA Kit

REF

DEIA050J



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The assay is intended for the quantitative determination of hemoglobin/haptoglobin complex in stool.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-Assay Precision: 2.49%-3.89%

Inter-Assay Precision: 2.64%-4.29%

Sensitivity

0.054 µg/mL