



User's Manual

HSV Type 1 rec. gG1 IgG- ELISA Kit



DEIA05535



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The HSV 1 recombinant gG1 IgG-ELISA allows the detection of HSV 1 infection in presence of antibodies to HSV 2 in human serum or plasma (citrate).

Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

The diagnostic sensitivity is >95 %.