



User's Manual

SPERM-ANTIBODY ELISA Kit



DEIA05703



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Sperm Antibody ELISA is a reliable and quantitative test for the determination antibodies directed against human spermatozoa. This test is intended for the use with serum.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

31-250 U/mL