



User's Manual

Toxocara ELISA Kit



DEIA05740



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the qualitative screening of serum IgG antibodies to Toxocara using an Enzyme-Linked Immunosorbent Assay (ELISA) technique.

Storage

Reagents, strips and bottled components:

Store at 2-8°C.

Squeeze bottle containing diluted wash buffer may be stored at room temperature. For more detailed information, please download the following document on our website.

Sensitivity

Sensitivity of 87.5% (21/24)