



User's Manual

B2GP1 ELISA Kit



DEIA05760



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

An enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of IgG or IgM antibodies to Beta 2 Glycoprotein 1 (B2GP1) in human serum as an aid in diagnosis of autoimmune thrombosis in patients with Systemic Lupus Erythematosus (SLE) or lupus like disorders in conjunction with other laboratory and clinical findings.

Storage

Store all reagents at 2-8°C. Do not freeze. For more detailed information, please download the following document on our website.

Detection Limit

The limits of detection (LoD) for these assays were determined to be 2.2 EU/mL for IgG and 5.3 EU/mL for IgM based on 60 replicates of the blank and 10 replicates each of 6 low-level(NHS) samples.

