



User's Manual

Hepatitis B Virus Nucleic Acid Related Antigen, HBV-NRAg ELISA Kit

REF

DEIA061



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This kit is an enzyme -linked immunosorbent assay (ELISA) for qualitative detection of HBV pre-S1 antigen and core antigen in hepatitis B virus nucleic acid related antigen (HBV NRAg) in human serum or plasma. It is intended for use in clinical laboratories for diagnosis and management of patients related to infection with hepatitis B virus.

Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

Human Anti-HBV Nucleic Acid Related Antigen ELISA Kit