



## User's Manual

# S100A12 ELISA Kit



DEIA062J



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The described Enzyme-Linked-Immuno-Sorbent-assay (ELISA) Kit is intended for the quantitative determination of S100A12 in stool.

### Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

### Precision

Intra-Assay Precision: 3.5%-7.1%

Inter-Assay Precision: 15.42%-20.18%

### Sensitivity

0.016 ng/mL