



User's Manual

Antibody to Treponema Pallidum (Syphilis) ELISA Kit



DEIA069



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Anti-TP ELISA is an enzyme-linked immunosorbent assay for qualitative determination of antibodies to Treponema Pallidum in human serum or plasma samples. The assay is intended to be used in clinical laboratories for diagnosis and management of patients related to infection with Treponema Pallidum.

Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.