



**User's Manual**

# Total Antibody to Hepatitis E Virus, HEV-Ab ELISA Kit



DEIA076



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This HEV-Ab ELISA kit is an enzyme linked-immunosorbent assay for in vitro qualitative detection of total antibodies (IgG, IgM, etc.) to hepatitis E virus in human or animal serum or plasma. It is intended to be used as an aid in supplementary diagnosis to hepatitis E infection, prevalence studies among populations as well as in zoonosis related research on hepatitis E virus.

### Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.

### Sensitivity

Human Anti-Hepatitis E Virus Total Antibody ELISA Kit