



## User's Manual

# Adiponectin Total ELISA Kit



DEIA078J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The described Enzyme-Linked-Immuno-Sorbent-assay (ELISA) Kit is intended for the quantitative determination of human adiponectin in serum and plasma.

### Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

### Precision

Intra-assay Precision: 3.2%-3.4%

Inter-assay Precision: 6.1%-6.3%

### Sensitivity

0.338 ng/mL