



User's Manual

Alpha-Fetoprotein, AFP ELISA Kit (Quantitative)



DEIA080



96T




This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This AFP ELISA kit is an enzyme linked immunosorbent assay (ELISA) for in vitro quantitative determination of α -fetoprotein (AFP) concentrations in the range of 2-400ng/mL in human serum or plasma samples.

Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.