



User's Manual

Carbonylated Proteins ELISA Kit



DEIA083J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Kit is intended for the determination of protein carbonyls in biological samples such as EDTA-plasma, bronchoalveolar lavage fluid and cerebrospinal fluid, cell extracts and other soluble protein samples.

Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

Precision

Intra-assay Precision: 5.80%-9.86%

Inter-assay Precision: 6.36%-9.72%

Sensitivity

20 pmol/mg