



User's Manual

Retinol binding Protein (RBP4) (1-point-calibration) ELISA Kit



DEIA091J



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA



Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com



Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Immundiagnostik assay is intended for the quantitative determination of free Retinol-binding protein (RBP)/RBP4 as well as RBP4 complexed with transthyretin in plasma, serum and urine.

Storage

Unopened Kit: Store at 2-8°C. Do not use past kit expiration date.

Precision

Intra-Assay Precision: 5%-5%

Inter-Assay Precision: 9.7%-9.8%

Sensitivity

0.9 µg/L